

ATLANTIC/ PACIFIC

Outdoor/indoor
power wheelchairs

*Heavy weight use
Perfect driving
characteristics
Maximum comfort
Fully range of adjust-
ment options*



Pacific / Atlantic. When extra weight capacity is required. Heavy weight makes special demands of a powered wheelchair. For this reason, Handicare has designed a modern electric wheelchair for people with a weight up to 200 kg as standard. Available as RWD (Pacific) and FWD (Atlantic), this wheelchair offers excellent driving characteristics for indoor and outdoor use. Uneven road surfaces will pass you by unnoticed due to the efficient suspension system. Halogen lights and indicators ensure safe driving at night, as you are clearly visible in the dark. Combining a proven concept with durable components and powerful motors specially developed for heavy weight use, the Pacific and Atlantic guarantee high reliability and durability. The functional and comfortable seating system is fully adjustable and offers various cushion types. A full range of electrical adjustment options is available, and of course, individual adjustments can be made to meet specific needs.

ATLANTIC/ PACIFIC



Heavy weight use

The Pacific and Atlantic are specially designed for heavier people. Reliability and durability are guaranteed.

Perfect driving characteristics

Highly manoeuvrable, the Pacific and Atlantic ensure a smooth ride, even on rugged street terrain. Powerful motors and excellent suspension. The sophisticated electronics systems can be adjusted to meet personal driving preferences.

Maximum comfort

The fully adjustable seating system provides maximum safety and comfort. Various cushions (width 60 or 65 cm) and leg rests, incl. length compensating leg rests, are available.

Fully range of adjustment options

The Pacific/ Atlantic can be fully customized to satisfy the specific requirements of the individual. A wide range of electrical adjustment options, incl. backrest recline, tilt-in-space and leg rests, is available.

DIMENSIONS			Atlantic		Pacific	
Axle configuration			FWD	FWD	RWD	RWD
Maximum user weight	kg	stone	200	32	200	32
Batteries	Ah	Ah	74	74	74	74
Total weight excl. batteries	kg	stone	85	13.4	85	13.4
Total weight incl. batteries	kg	stone	115	18	115	18
Total length	mm	inch	1105	43.5	1105	43.5
Total width	mm	inch	740	29.1	740	29.1
Seat angle	°	°	0 - 10	0 - 10	0 - 10	0 - 10
Seat depth	mm	inch	415 - 540	16.3 - 21.3	415 - 540	16.3 - 21.3
Seat width	mm	inch	600 - 730	23.6 - 28.7	600 - 730	23.6 - 28.7
Seat height exc. Cushion	mm	inch	440 - 465	17.3 - 18.3	440 - 465	17.3 - 18.3
Seat height incl. cushion	mm	inch	510 - 535	20.1 - 21.1	510 - 535	20.1 - 21.1
Backrest angle	°	°	90 - 115	90 - 115	90 - 115	90 - 115
Back height	mm	inch	540	21.3	540	21.3
Diameter castor wheels (indoor)	mm	inch	250	9.8	220	8.7
Diameter castor wheels (outdoor)	mm	inch	250	9.8	220	8.7
Diameter drive wheels	mm	inch	350	13.75	350	13.75
PERFORMANCE						
Maximum speed	km/h	mph	8	5	8	5
Range (ISO 7176-4)	km	miles	25	15.5	25	15.5
Obstacle height	mm	inch	60	2.4	60	2.4
Obstacle height with kerb climber	mm	inch	-	-	-	-
Ground clearance	mm	inch	80	3.2	80	3.2
Turning radius (ISO 7176-5)	mm	inch	700	27.6	990	39
ELECTRICAL OPTIONS						
Electrical backrest adjustment	°	°	-	-	-	-
Electrical seat tilt adjustment	°	°	-	-	-	-
Electrical high low	mm	inch	-	-	-	-
Electrical leg rests	°	°	-	-	-	-

EN12182

EN12184



handicare

The Netherlands
T +31 (0)492 593 888
F +31 (0)492 537 931
info@handicare.com
www.handicare.com

UK
T +44 (0)1932 858 687
F +44 (0)1932 848 384
info@handicare.co.uk
www.handicare.co.uk