



Yes, you can.[®]

 **Invacare[®] Leo[™]**

Safe, solid, secure: the *Invacare[®] Leo[™]*.

The *Invacare[®] Leo[™]* is a 4-wheel scooter designed for all those that value their independence and wish to get out and about unaided. Safety is a key feature of the **Leo[™]**, but this does not detract from its stylish and sporty looks. **Leo[™]** offers users the freedom and confidence to enjoy their essential daily outings and leisure excursions. With exactly high standards of design and safety engineering, the **Leo[™]** undoubtedly sets a new market standard.





➤ **For a safe, powerful drive**

The *Invacare*® **Leo**™ offers a smooth, safe drive in a scooter that is easy to control and manoeuvre. The **Leo**™ intelligent lighting package includes indicators that automatically switch off and a rear brake light for safe driving. The lighting system is designed so that if a bulb should blow, all other lights still continue to function. In addition to a visual warning for the battery level display, users are also alerted with an audible warning (beep) when the battery needs recharging.

➤ **Compact**

The **Leo**™ is a compact mini-scooter that can be easily dismantled for transportation during longer trips or weekend breaks. With a smooth drive and easy manoeuvrability, getting round tight corners in the local shopping centre, museum or leisure complex becomes refreshingly easy.

➤ **Reliable**

Users can trust the **Leo**™ to take them wherever they need to go. From a trip to the shops or a ride to the park, the **Leo**™ is rendered sturdy and reliable thanks to its exceptionally high quality electronic engineering and power. These elements are protected from corrosion and water, (splash protection) which safeguard the vital workings of the scooter against dirt and water. A two-step disengaging lever protects the scooter from free-wheeling if inadvertently knocked out of drive mode.

➤ **Easy access**

Access to the battery and electronics for maintenance or servicing is much easier with the **Leo**™ as only one part needs to be removed. This simple system also facilitates removal of batteries when disassembling the **Leo**™ for transport.



Features



Finite tiller movement

Users can easily adjust the tiller to suit their needs.




Brake light

Rear brake lighting warns other people that the scooter is slowing down, thus preventing collisions.



Comfort seating

Swivel functionality and sliding rails ensure optimal seating position.



Two-step disengaging lever

This unique concept prevents the scooter from accidentally freewheeling.



Built in anti-splash guards

Protects electronics and the transaxle from dirt and water.



Disassembly into five parts

No tools are required to dismantle
Invacare® Leo™

Options & Accessories



Crutch & cane holder

Keeps your crutch or cane stable during the drive.



Lockable front box

Secure box for personal valuables.



Seatbelt

Keeps the user well-positioned in the seat.



Large mirror

Can be mounted left or right.



Heavy-duty storage cover

Keeps the scooter clean and dry.



Charging off-board

Simply remove batteries to charge with the off-board charger kit.

Technical Data



47 cm




41 cm



46 - 56 cm




47,5 cm



59 cm



122 cm




59 kg (excl. batteries)
83 kg (incl. batteries)



136 kg



2 x 36 Ah



131 cm



6 cm



7 cm



10°



240 W




36 km



8 km/h

Colour of shrouds




Moonstone Silver



Sapphire Blue

Upholstery



Grey/Black vinyl



All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.



INVACARE® Limited
South Road · Bridgend Industrial Estate · BRIDGEND CF31 3PY
Tel. +44 (0) 1656 647327 · Fax +44 (0) 1656 649016
E-mail: uk@invacare.com · www.invacare.co.uk
Sales Order E-mail: ordersuk@invacare.com

INVACARE® Ireland Ltd
Unit 5 · Seatown Business Campus
Seatown Road · SWORDS · County Dublin
Tel. (+353) 181 107 84 · Fax (+353) 181 107 85
Email: ireland@invacare.com · www.invacare.ie