



Yes, you can.[®]

Invacare[®] Alu Lite

A lighter way to get around

The *Invacare[®] Alu Lite* is the basic-entry wheelchair in the Invacare[®] manual range. The lightweight aluminium construction of this transit wheelchair makes it easy to manoeuvre, even for an elderly carer, as well as providing a wheelchair that is easily transportable.



NEW

➤ Easy to transport and store

Being lightweight (aluminium chassis), compact and with a foldable back, the **Invacare® Alu Lite** is very easy to handle, transport and store. Simply remove the footrests, fold the back and chassis and store away.



➤ Pushing? No problem!

The overall manoeuvrability of a wheelchair is crucial - which is why the **Alu Lite** is designed to be very easy to push and manoeuvre, both indoors and out. Due to the smaller overall dimensions and general compactness of this model, getting around in tight spaces is never a problem.

Options & Features



Folding backrest

Reduces overall size when folded for transportation or storage.



Swing away footrest

Footplates fold up as shown, and footrests swing away, allowing easier access in or out of the chair.



Foldable cross frame

Lightweight frame folds flat for easier portability.

Upholstery



Black nylon

Frame colour



Blue

Technical data

40,5 / 45,5 cm	43 cm	49 cm	45 cm	23 cm
34 - 46 cm	3°	0°	52 / 57 cm	94,5 cm
103 cm	76 cm	12 kg	max. 100 kg	10,2 kg

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.



INVACARE® Limited
 South Road · Bridgend Industrial Estate · BRIDGEND CF31 3PY
 Tel. +44 (0) 1656 647327 · Fax +44 (0) 1656 649016
 E-mail: uk@invacare.com · www.invacare.co.uk
 Sales Order E-mail: ordersuk@invacare.com

INVACARE® Ireland Ltd
 Unit 5 · Seatown Business Campus
 Seatown Road · SWORDS · County Dublin
 Tel. (353) 181 107 84 · Fax (353) 181 107 85
 Email: ireland@invacare.com · www.invacare.ie