



Yes, you can.®

 Invacare® **XLT**

### Get Perfect Support - Remain Active!

The new Invacare® **XLT** is the ideal active, rigid chair for day to day use. Designed to be perfectly adjusted to fit the body the **XLT** offers its users a wide range of various levels of support, enabling and even improving independence and mobility.

**NEW**





### ➤ **Simply Compatible**

Complex chair-specific solutions and spare parts are a thing of the past with the new *Invacare*® **XLT**. Benefit from the extensive and readily available range of *Invacare*® options which can be fitted to the **XLT** without any hassle.

### ➤ **Stability, Rigidity and Security for an Active Life**

Made entirely from titanium the **XLT** gives its user stability and rigidity without compromising on weight. In combination with the wide range of options, features and configuration possibilities offered with the **XLT** the user will benefit from a chair which offers the maximum security and support needed when going about busy day to day life.

### ➤ **The Perfect Fit**

The **XLT** has been designed with the individual in mind. Measured, constructed and fitted to suit the precise requirements of each user, the **XLT** is the ideal chair to help regain, retain and improve both mobility and independence.

### ➤ **Proven Quality**

*Invacare*®, together with *Rea*™ have made no compromises when it comes to the new **XLT**. Designed and built at the *Rea* factory in Sweden, users will benefit from the strong tradition of both *Invacare*® and *Rea*™ quality fused to create one unique chair: The *Invacare*® **XLT**.



## Options



### Swing Away

for easy transfer and maneuverability as well as comfort.



**Front Castor Attachments**  
for stability and performance.



### Titanium Foottube

to improve comfort as well as lightness.



### Tyres

① solid, ② pneumatic profiled or puncture proof, ③ pneumatic non marking, ④ low profile, ⑤ pneumatic profiled or puncture proof, ⑥ black marathon evolution.



### Frontwheels

① 75 x 24 hard (3"), ② 100 x 27 medium hard (4"), ③ 125 x 27 medium hard (5"), ④ 140 x 37 soft (5.5"), ⑤ 150 x 27 medium hard (6"), ⑥ 150 x 30 soft (6"), ⑦ 150 x 32 pneumatic (6"), ⑧ 180 x 40 pneumatic (7"), ⑨ 180 x 40 soft (7")

## Options



### Armrests

Wide range of armrests for best support.



### Anti tipper device


The Anti-tip device can easily be folded away by the foot. The device is also angled outwards, to give space for the feet for the attendant.



### Sporty Wheel

with black Marathon Evolution Plus tyres for design.

## Technical data




35,5, 38, 40,5, 43, 45,5, 48cm / 50,5 cm Swing



40 / 45 / 50 (on 4 cm adj.)



40 - 52 cm



20 - 49 cm




36 - 48 cm  
27 - 35 cm high mounted



0° - 6°



-10° - +10°




Seat width + 20 cm



76 - 113 cm



Active 83 - 103 cm,  
Dynamic 73 - 88 cm,  
Swing 82 - 123 cm



9,8 kg



max 135 kg




Transport weight:  
6 kg



Crash-tested  
ISO 7176-19\*

\* The Invacare® **XLT** complies with ISO 7176-19 specifications. Note that special conditions apply regarding accessories. For detailed information regarding vehicular transportation for this wheelchair, please consult the following pages on our website: [www.invacare-rea.com](http://www.invacare-rea.com) under Rea wheelchairs/Technical downloads.

## Frame colours



Blue



Red



Black



Grey



Blasted Titanium\*  
not for Swing



Orange\*

\* Available as special orders

## Upholstery



Black Nylon



Black Jemima

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.



INVACARE® Rea AB  
Box 200 · Växjövägen 303  
SE-343 75 DIÖ · Sweden  
Tel: +46 476 535 00 · Fax + 46 476 535 99  
[www.invacare-rea.com](http://www.invacare-rea.com)

