

**QUANTUM**  
**MOVING FORWARD**

A Division of Pride Mobility Products® Ltd.



# Artemis

Rear-Wheel Drive Power Chair

Artemis shown with Vogue seating  
and PG VR2 controller

## The Artemis Rear-Wheel Drive Power Chair

The Artemis is designed for clients who demand full rehab capability and tight space manoeuvrability without sacrificing traditional rear-wheel drive performance.



# Specifications

<b>DRIVE WHEELS</b>	30.5 cm (12")
<b>CASTORS</b>	20 cm (8") solid
<b>ANTI-TIPS</b>	5" (13 cm) rear
<b>MAXIMUM SPEED</b>	Up to 6 km/h (4 mph) <sup>1</sup>
<b>BRAKES</b>	Intelligent Braking (electronic, regenerative disc brakes)
<b>GROUND CLEARANCE</b>	7.6 cm (3")
<b>OVERALL SIZE<sup>2</sup></b>	Length: 78 cm (30.7") Width: 61 cm (24")
<b>SEATING</b>	Sling-back with solid baseboard / Comfort seating
<b>DRIVE TRAIN</b>	2-motor, rear-wheel drive
<b>BATTERY REQUIREMENTS</b>	Two 12 volt, deep cycle
<b>BATTERY SIZE</b>	50 AH (standard) 60 AH (optional)
<b>BATTERY CHARGER</b>	Off-board, 8 amp (standard)
<b>AVAILABLE ELECTRONICS</b>	60 amp, PG VR2
<b>WEIGHT CAPACITY</b>	136 kg (21 Stone) maximum (user weight)
<b>BATTERY WEIGHT<sup>3</sup></b>	<b>50 AH:</b> 17 kg (38 lbs.) (each) <b>60 AH:</b> 19.5 kg (43 lbs.) (each)
<b>WARRANTY</b>	Lifetime limited warranty on frame 1-year limited warranty on electronics 18-month warranty on drive control system 2-year limited warranty on drive train

# Colours



Metallic Red    Metallic Blue    Metallic Silver    Metallic Black    Sunburst Orange

The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.

# Features

- Outstanding tight-space manoeuvrability
- Narrow 61 cm (24") width and compact 78 cm (30.7") length
- Maximum speed of up to 6 km/h (4 mph)<sup>1</sup>
- Adjustable seat frame
- Unique new design, seating system is integrated with the base
- Width and depth adjustable armrests
- Glide Motors for added torque

# Options

- Flat-free tyres
- Power tilt
- Power recline
- Power tilt and recline
- Pediatric legrests
- TRU-Comfort® cushions
- Lap belt
- 70° swing-away footrests
- 80° heavy-duty leg rests
- 70° elevating leg rests
- Power elevating leg rests
- Heel loops
- LED lighting

### FOOTNOTES:

- 1) Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3 %.
- 3) Battery weight may vary based on manufacturer.



A Division of Pride Mobility Products® Ltd.

32 Wedgwood Road, Bicester,  
Oxfordshire OX26 4UL  
(Tel) 01869 324 600  
(Fax) 01869 323 070  
www.quantumrehab.co.uk

