

# POWERTEC F16



## Special Features

- The all-in-one power unit can be easily mounted and dismantled to the manual chair giving you the option of a manual or power chair as and when you want it
- Compatible with the Quickie 2, Quickie 2 Millennium and RXS manual wheelchair
- Adaptable to a number of chair widths - 14" to 18" frames
- Comes with 25AH maintenance free batteries and a 6 amp battery charger
- Power unit and chair can be easily dismantled and transported
- Appropriate for both indoor and outdoor use



transforms from  
power to manual  
in seconds



**SUNRISE**

*Improving people's lives*

## Features and Accessories

Maintenance free batteries	●	Automatic charger	●
Swing-away legrests	●	Height adjustable armrests	●
Puncture proof castor wheels	○	Pneumatic castor wheels	●
Rear wheel drive	●	Sling upholstery	●
Padded upholstery	○	Remote pod	●
Folding frame	●	Attendant Control	○
Dual control	○	Kerb climber	○
Puncture proof tyres	○	Lights and indicators	○
Headrest	○	Elevated legrest	○
Lap strap	○	Calf strap	○
Manual handbrakes	○	Thoracic supports	○
Foam cushions	○	Toe straps	○
Weatherproof cape	○	Angle adjustable footplates	○

Sunrise reserves the right to change specifications without prior notice.

● Standard ○ Optional

## Performance, weights and measures

Maximum range	13km
Kerb climbing*	10cm
Max. forward speed	6kph
Max. user weight	118kg
Max. gradient	25%
Overall length	115cm
Overall width	64cm
Folded width	33cm
Total weight (inc. standard seating and batteries)	53kg
Weight of heaviest part	19kg
Seat depth	36-46cm
Seat width	36-46cm
Seat height	48cm
Backrest height	43cm
Backrest angle	8°

Temperature, user weight, terrain and battery capacity/age may effect overall performance.

Weights and measures are approximate and should only be used as a guide.

\* When optional kerb climber is fitted

## Colour options



Two powerpacks fit perfectly into place in the power module



Height adjustable armrests to suit the user



Compact power module can be easily dismantled



Swing-away legrests

CE All products are CE marked to the MDD 93/42/EEC - 14/06/93. Radio interference is under EEC Directive 89/336/EEC and battery chargers are under EEC Directives 73/23/EEC and 89/336/EEC.



**SUNRISE**

Sunrise Medical Limited

High Street Wollaston West Midlands DY8 4PS England

Tel: 01384 44 66 88 Fax: 01384 44 66 99