

QUICKIE®

Wide range of sizes and options that meet clinical needs...

FLYQR



- > VR2 control system
- > Wide range of seat sizes, from 13" wide - ideal for paediatrics and juniors
- > Suitable for Indoor + Outdoor use



Technical details:



The control system can easily be programmed to suit your individual performance requirements



The ergonomically designed "Y" shaped backrest ensures total comfort







New height adjustable detachable armrests



Twin 30Ah (std) or 40Ah batteries



The F4oVR disassembles without tools for easy transport



Max. range: ISO 7176 - 4

23 km

Turning radius:

1680 mm

Speed:

6 km/h

Battery:

Twin 30Ah or 40Ah

Kerb climbing

max. 5 cm / max. 10 cm with kerb climber

Max. safe slope

14°

* TRANSIT: The Quickie F4oVR is suitable for user transportation as a seat in a vehicle as long as it is used in conjunction with the standard approved tie down system and occupant restraint system.

STEPS
Sunrise Training and Education Programmes

SUNRISE MEDICAL

Sunrise Medical Limited
High Street
Wollaston, Stourbridge
West Midlands DY8 4PS, England
Phone +44 (0) 1384 446666
Fax +44 (0) 1384 446644
www.sunrisemedical.com