


# Science of Seating

JAY has over thirty years of experience and research in clinical seating including client research, anthropometric databases, pressure mapping and material science.



Jay Lite has been developed using the latest materials, profiling and clinical research in providing the ideal cushion for the active user.

**Microclimatic Material**



**Tension Reduction**



## Performance and Easy transfers

Weighing from just 780 grams at **40 x 40 cm** cushion the **JAY LITE** is **one of the lightest** pressure re-distribution cushions available.

**Reaction Force Reduction**



**Pre-contouring**



**Firmness Layering**



**Positioning Capability**



JAY LITE also covers the other key issues for the active client.

## Superior Pressure Distribution

JAY LITE uses ischial scoop cut-outs combined with the anthropometrically researched Optiwell™ technology to redirect the pressure away from the bony prominences of the pelvis and onto the hips and legs.

Firm Oletex® and soft EVA foam layering ensures comfort without sacrificing stability, essential for the active client.

## Positioned for Action

JAY LITE comes in 44 sizes and two versions. The standard base is ideal for comfort, pressure redistribution, heat and moisture dissipation for most clients. For those that require more stability the positioning option provides lateral and medial thigh supports.

## Superior Heat and Moisture Dissipation

An active life creates heat generation and a need to stay “cool”. JAY LITE uses the new 3DX™ microclimatic “over-size” covers and an Olotex® “airflow” base with air vents and porous foams.

## Stay Cool to Look “Cool”

JAY LITE not only ensures you are comfortable and well poised for your day, it looks good too. Key features including the non-slip base, backrest notches and easy-grip transfer handle to ensure the JAY LITE stays in place and is easy to transfer.

## Technical Specifications

44 seat sizes from 35 x 35 cm to 60 x 60 cm

Maximum User weight: 113 kg

Cushion weight 780 grams for 40 x 40 cm standard profile

Cushion height: 7 cm



We would like to thank Tarek Rasouli for his participation in this brochure.

Sunrise Medical Limited  
High Street  
Wollaston, Stourbridge  
West Midlands, DY8 4PS  
T +44 (0) 1384 446666  
F +44 (0) 1384 446644  
[www.sunrisemedical.com](http://www.sunrisemedical.com)



**JAY<sup>®</sup> LITE** – feel the LITE

**LIGHT** from only 780 grams\* with Olotex® base  
**COOL** moisture and heat dissipation with 3DX™ spacer fabric  
**ACTIVE** 'positioned to go' with Optiwell™ technology  
plus 44 sizes and positioning options

Want more information or details of your nearest dealer?  
... register now on [www.jay-seating.com](http://www.jay-seating.com)

\*cushion dimension 40x40 cm, and cushion 35x35 cm is even lighter



Reaction Force Reduction




Pre-contouring



Firmness Layering



Microclimatic Material



Positioning Capability



Tension Reduction





**JAY<sup>®</sup> LITE** – stay cool